

Product datasheet for **GA108267**

Human SPECCIL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPECCIL CRISPRa kit – CRISPR gene activation of human sperm antigen with calponin homology and coiled-coil domains 1 like
Symbol:	SPECCIL
Components:	GA108267G1, SPECCIL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA108267G2, SPECCIL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA108267G3, SPECCIL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q69YQ0</u>
Synonyms:	CYTSA; GBBB2; OBLFC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

