

## Product datasheet for **GA108270**

### Human SIK3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | SIK3 CRISPRa kit - CRISPR gene activation of human SIK family kinase 3  |
| Symbol:       | SIK3  |
| UniProt ID:   | <a href="#">Q9Y2K2</a>  |
| Synonyms:     | L19; QSK; SIK-3   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)