

Product datasheet for **GA108276**

Human PIP5KI gamma (PIP5K1C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIP5K1C CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4-phosphate 5-kinase type 1 gamma
Symbol:	PIP5K1C
UniProt ID:	O60331
Synonyms:	LCCS3; PIP5K-GAMMA; PIP5K1-gamma; PIP5Kgamma
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)