

Product datasheet for **GA108278**

Human PPWD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPWD1 CRISPRa kit - CRISPR gene activation of human peptidylprolyl isomerase domain and WD repeat containing 1
Symbol:	PPWD1
UniProt ID:	<u>Q96BP3</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)