

## Product datasheet for **GA108293**

### Human **KCNH4** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNH4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 4
Symbol:	KCNH4
UniProt ID:	<a href="#">Q9UQ05</a>
Synonyms:	BEC2; ELK1; Kv12.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)