

Product datasheet for **GA108294**

Human **KCNH3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNH3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 3
Symbol:	KCNH3
UniProt ID:	<u>Q9ULD8</u>
Synonyms:	BEC1; ELK2; Kv12.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)