

Product datasheet for **GA108298**

Human **ITGB3BP** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	ITGB3BP CRISPRa kit - CRISPR gene activation of human integrin subunit beta 3 binding protein
Symbol:	ITGB3BP
UniProt ID:	Q13352
Synonyms:	CENP-R; CENPR; HSU37139; NRIF3; TAP20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)