

## Product datasheet for **GA108299**

### Human TMED3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | TMED3 CRISPRa kit - CRISPR gene activation of human transmembrane p24 trafficking protein 3 |
| Symbol:       | TMED3   |
| UniProt ID:   | <a href="#">Q9Y3Q3</a>  |
| Synonyms:     | C15orf22; P24B; p26   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                     |

#### Product images:



[View online >](#)