

Product datasheet for **GA108313**

Human **SLC35A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member A3
Symbol:	SLC35A3
UniProt ID:	Q9Y2D2
Synonyms:	AMRS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)