

Product datasheet for **GA108332**

Human QPRT activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | QPRT CRISPRa kit - CRISPR gene activation of human quinolinate phosphoribosyltransferase |
| Symbol: | QPRT |
| UniProt ID: | Q15274 |
| Synonyms: | HEL-S-90n; QPRTase |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)