

Product datasheet for **GA108338**

Human TGDS activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TGDS CRISPRa kit – CRISPR gene activation of human TDP-glucose 4,6-dehydratase
Symbol:	TGDS
Components:	GA108338G1, TGDS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108338G2, TGDS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108338G3, TGDS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O95455</u>
Synonyms:	CATMANS; SDR2E1; TDPGD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

