


Product datasheet for **GA108351**

Human **LRRC8B** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC8B CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 8 VRAC subunit B
Symbol:	LRRC8B
Components:	GA108351G1, LRRC8B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108351G2, LRRC8B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108351G3, LRRC8B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P9F7</u>
Synonyms:	TA-LRRP; TALRRP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)