

## Product datasheet for **GA108362**

### Human R3HDM1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | R3HDM1 CRISPRa kit - CRISPR gene activation of human R3H domain containing 1                                                                                                                                                                                                                                  |
| Symbol:       | R3HDM1                                                                                                                                                                                                                                                                                                        |
| Components:   | GA108362G1, R3HDM1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108362G2, R3HDM1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108362G3, R3HDM1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q15032</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | R3HDM                                                                                                                                                                                                                                                                                                         |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)