

## Product datasheet for **GA108368**

### Human SRRM2 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SRRM2 CRISPRa kit - CRISPR gene activation of human serine/arginine repetitive matrix 2
Symbol:	SRRM2
UniProt ID:	<a href="#">Q9UQ35</a>
Synonyms:	300-KD; Cwc21; CWF21; HSPC075; SRL300; SRm300
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)