

Product datasheet for **GA108376**

Human PIK3R5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3R5 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase regulatory subunit 5
Symbol:	PIK3R5
UniProt ID:	<u>Q8WYR1</u>
Synonyms:	F730038I15Rik; FOAP-2; p101; P101-PI3K
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)