


Product datasheet for **GA108377**

Human TNPO3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNPO3 CRISPRa kit - CRISPR gene activation of human transportin 3
Symbol:	TNPO3
Components:	GA108377G1, TNPO3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108377G2, TNPO3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108377G3, TNPO3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y5L0</u>
Synonyms:	IPO12; LGMD1F; MTR10A; TRN-SR; TRN-SR2; TRNSR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)