

## Product datasheet for **GA108380**

### Human SLC16A8 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC16A8 CRISPRa kit - CRISPR gene activation of human solute carrier family 16 member 8
Symbol:	SLC16A8
Components:	GA108380G1, SLC16A8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108380G2, SLC16A8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108380G3, SLC16A8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">O95907</a></u>
Synonyms:	MCT3; REMP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)