

Product datasheet for **GA108388**

Human **TTC33** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC33 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 33
Symbol:	TTC33
UniProt ID:	Q6PID6
Synonyms:	OSRF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)