

Product datasheet for **GA108390**

Human PSD4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSD4 CRISPRa kit - CRISPR gene activation of human pleckstrin and Sec7 domain containing 4
Symbol:	PSD4
UniProt ID:	Q8NDX1
Synonyms:	EFA6B; TIC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)