

Product datasheet for **GA108434**

Human PHLDA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHLDA3 CRISPRa kit - CRISPR gene activation of human pleckstrin homology like domain family A member 3
Symbol:	PHLDA3
UniProt ID:	Q9Y5J5
Synonyms:	TIH1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)