

Product datasheet for **GA108436**

Human SH3BP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3BP1 CRISPRa kit - CRISPR gene activation of human SH3 domain binding protein 1
Symbol:	SH3BP1
UniProt ID:	Q9Y3L3
Synonyms:	ARHGAP43
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)