

Product datasheet for **GA108437**

Human TSSK2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSSK2 CRISPRa kit - CRISPR gene activation of human testis specific serine kinase 2
Symbol:	TSSK2
UniProt ID:	Q96PF2
Synonyms:	DGS-G; SPOGA2; STK22B; TSK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)