

Product datasheet for **GA108446**

Human KCNE1L (KCNE5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNE5 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily E regulatory subunit 5
Symbol:	KCNE5
UniProt ID:	Q9UJ90
Synonyms:	KCNE1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)