


Product datasheet for **GA108458**

Human PLD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLD3 CRISPRa kit - CRISPR gene activation of human phospholipase D family member 3
Symbol:	PLD3
UniProt ID:	Q8IV08
Synonyms:	AD19; HU-K4; HUK4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)