

Product datasheet for **GA108459**

Human ARFIP2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARFIP2 CRISPRa kit - CRISPR gene activation of human ADP ribosylation factor interacting protein 2
Symbol:	ARFIP2
Components:	GA108459G1, ARFIP2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108459G2, ARFIP2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108459G3, ARFIP2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P53365</u>
Synonyms:	POR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)