

Product datasheet for **GA108467**

Human ZKSCAN5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZKSCAN5 CRISPRa kit - CRISPR gene activation of human zinc finger with KRAB and SCAN domains 5
Symbol:	ZKSCAN5
Components:	GA108467G1, ZKSCAN5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108467G2, ZKSCAN5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108467G3, ZKSCAN5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y2L8</u>
Synonyms:	ZFP95; ZNF914; ZSCAN37
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)