


Product datasheet for **GA108472**

Human SH3BP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3BP4 CRISPRa kit - CRISPR gene activation of human SH3 domain binding protein 4
Symbol:	SH3BP4
Components:	GA108472G1, SH3BP4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108472G2, SH3BP4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108472G3, SH3BP4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9P0V3</u>
Synonyms:	BOG25; TTP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)