

Product datasheet for **GA108473**

Human **SGK3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SGK3 CRISPRa kit - CRISPR gene activation of human serum/glucocorticoid regulated kinase family member 3
Symbol:	SGK3
UniProt ID:	<u>Q96BR1</u>
Synonyms:	CISK; SGK2; SGKL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)