

Product datasheet for **GA108475**

Human PKC nu (PRKD3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKD3 CRISPRa kit - CRISPR gene activation of human protein kinase D3
Symbol:	PRKD3
UniProt ID:	O94806
Synonyms:	EPK2; nPKC-NU; PKC-NU; PKD3; PRKCN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)