

## Product datasheet for **GA108477**

### Human KCNE4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNE4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily E regulatory subunit 4
Symbol:	KCNE4
UniProt ID:	<a href="#">Q8WWG9</a>
Synonyms:	MIRP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)