

## Product datasheet for **GA108484**

### Human EID1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	EID1 CRISPRa kit - CRISPR gene activation of human EP300 interacting inhibitor of differentiation 1
Symbol:	EID1
Components:	GA108484G1, EID1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108484G2, EID1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108484G3, EID1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q9Y6B2</a></u>
Synonyms:	C15orf3; CRI1; EID-1; IRO45620; PNAS-22; PTD014; RBP21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)