

Product datasheet for **GA108488**

Human SDF2L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SDF2L1 CRISPRa kit - CRISPR gene activation of human stromal cell derived factor 2 like 1
Symbol:	SDF2L1
Components:	GA108488G1, SDF2L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108488G2, SDF2L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108488G3, SDF2L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HCN8</u>
Synonyms:	AP000553.C22.4; dihydropyrimidinase-like 2; OTTHUMT00000075032; PWP1-interacting protein 8; SDF2-like 1; stromal cell-derived factor 2-like 1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)