

Product datasheet for **GA108494**

Human IL17RA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL17RA CRISPRa kit - CRISPR gene activation of human interleukin 17 receptor A
Symbol:	IL17RA
UniProt ID:	Q96F46
Synonyms:	CANDF5; CD217; CDw217; hIL-17R; IL-17RA; IL17R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)