

Product datasheet for **GA108500**

Human ARHGAP8 activation kit by CRISPRa

Product data:

| | |
|---------------|---------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | ARHGAP8 CRISPRa kit - CRISPR gene activation of human Rho GTPase activating protein 8 |
| Symbol: | ARHGAP8 |
| UniProt ID: | P85298 |
| Synonyms: | BPGAP1; PP610 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)