

## Product datasheet for **GA108507**

### Human IFIT5 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	IFIT5 CRISPRa kit - CRISPR gene activation of human interferon induced protein with tetratricopeptide repeats 5
Symbol:	IFIT5
UniProt ID:	<a href="#">Q13325</a>
Synonyms:	ISG58; P58; RI58
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)