

Product datasheet for **GA108518**

Human TP53TG3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TP53TG3 CRISPRa kit - CRISPR gene activation of human TP53 target 3
Symbol:	TP53TG3
UniProt ID:	<u>Q9ULZ0</u>
Synonyms:	P53TG3; TP53TG3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)