

Product datasheet for **GA108531**

Human RIPK5 (DSTYK) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | DSTYK CRISPRa kit - CRISPR gene activation of human dual serine/threonine and tyrosine protein kinase |
| Symbol: | DSTYK |
| UniProt ID: | Q6XUX3 |
| Synonyms: | CAKUT1; DustyPK; HDCMD38P; RIP5; RIPK5 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)