

Product datasheet for **GA108533**

Human RAB3GAP2 activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: RAB3GAP2 CRISPRa kit - CRISPR gene activation of human RAB3 GTPase activating non-catalytic protein subunit 2

Symbol: RAB3GAP2


Components: GA108533G1, RAB3GAP2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)
GA108533G2, RAB3GAP2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)
GA108533G3, RAB3GAP2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)
1 CRISPRa-Enhancer vector, SKU GE100056
1 CRISPRa scramble vector, SKU GE100077

UniProt ID: [Q9H2M9](#)

Synonyms: p150; RAB3-GAP150; RAB3GAP150; SPG69; WARBM2

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)