

Product datasheet for **GA108543**

Human BRI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	BRI3 CRISPRa kit - CRISPR gene activation of human brain protein I3
Symbol:	BRI3
Components:	GA108543G1, BRI3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108543G2, BRI3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108543G3, BRI3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O95415</u>
Synonyms:	I3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)