

Product datasheet for **GA108544**

Human ZNF324 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF324 CRISPRa kit - CRISPR gene activation of human zinc finger protein 324
Symbol:	ZNF324
Components:	GA108544G1, ZNF324 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108544G2, ZNF324 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108544G3, ZNF324 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O75467</u>
Synonyms:	ZF5128; ZNF324A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

