

Product datasheet for **GA108557**

Human TAF A5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF A5 CRISPRa kit - CRISPR gene activation of human TAF A chemokine like family member 5
Symbol:	TAF A5
UniProt ID:	Q7Z5A7
Synonyms:	QLLK5208; TAF A-5; TAF A5; UNQ5208
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)