

Product datasheet for **GA108580**

Human PREI3 (MOB4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MOB4 CRISPRa kit - CRISPR gene activation of human MOB family member 4, phocoin
Symbol:	MOB4
UniProt ID:	Q9Y3A3
Synonyms:	2C4D; CGI-95; MOB1; MOB3; MOBKL3; PHOCN; PREI3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)