

Product datasheet for **GA108585**

Human TECPR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TECPR1 CRISPRa kit - CRISPR gene activation of human tectonin beta-propeller repeat containing 1
Symbol:	TECPR1
Components:	GA108585G1, TECPR1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108585G2, TECPR1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108585G3, TECPR1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7Z6L1
Synonyms:	DKFZP434B0335
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)