

Product datasheet for **GA108599**

Human LETMD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LETMD1 CRISPRa kit - CRISPR gene activation of human LETM1 domain containing 1
Symbol:	LETMD1
Components:	GA108599G1, LETMD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108599G2, LETMD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108599G3, LETMD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P1Q0</u>
Synonyms:	1110019O13Rik; DKFZp586A011; HCCR-1; HCCR-2; HCCR1; HCRR-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)