

## Product datasheet for **GA108606**

### Human WDR51A (POC1A) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | POC1A CRISPRa kit - CRISPR gene activation of human POC1 centriolar protein A                                                                                                                                                                                                                                 |
| Symbol:       | POC1A                                                                                                                                                                                                                                                                                                         |
| Components:   | GA108606G1, WDR51A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108606G2, WDR51A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108606G3, WDR51A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8NBT0</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | PIX2; SOFT; WDR51A                                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)