

Product datasheet for **GA108607**

Human ZNF473 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF473 CRISPRa kit - CRISPR gene activation of human zinc finger protein 473
Symbol:	ZNF473
Components:	GA108607G1, ZNF473 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108607G2, ZNF473 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108607G3, ZNF473 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WTR7</u>
Synonyms:	ZFP100; ZN473
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)