

Product datasheet for **GA108611**

Human PLEKHG4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHG4 CRISPRa kit - CRISPR gene activation of human pleckstrin homology and RhoGEF domain containing G4
Symbol:	PLEKHG4
UniProt ID:	Q58EX7
Synonyms:	ARHGEF44; PRTPHN1; SCA4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)