

## Product datasheet for **GA108622**

### Human TMEM158 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM158 CRISPRa kit - CRISPR gene activation of human transmembrane protein 158 (gene/pseudogene)
Symbol:	TMEM158
UniProt ID:	<a href="#">Q8WZ71</a>
Synonyms:	BBP; p40BBP; RIS1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)