

## Product datasheet for **GA108625**

### Human C1orf43 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | C1orf43 CRISPRa kit - CRISPR gene activation of human chromosome 1 open reading frame 43 |
| Symbol:       | C1orf43  |
| UniProt ID:   | <a href="#">Q9BWL3</a>   |
| Synonyms:     | HSPC012; NICE-3; NS5ATP4; S863-3   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                  |

#### Product images:



[View online >](#)