

## Product datasheet for **GA108639**

### Human HDPTP (PTPN23) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PTPN23 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 23
Symbol:	PTPN23
UniProt ID:	<a href="#">Q9H3S7</a>
Synonyms:	HD-PTP; HDPTP; PTP-TD14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)